DNA Mixture Interpretation Principles: Observations from a NIST Scientific Foundation Review AAFS 2019 Workshop #10 (February 18, 2019; Baltimore, MD)



Mixture Interpretation Consistency Training Initiative

Dr. Ray Wickenheiser

New York State Police Crime Lab System Director






Outline

- Initiative to improve consistency of mixture interpretation
- Improvements made to mixture interpretation procedure lead by an internal mixture working group
- Training of new analysts coupled with retraining of existing experts
- Use in extensive examples in training and in manuals
- All analysts interpret the same complex 3 person mixtures
- Examine performance and seek opportunities for improvement

Mixture Training and Procedures

- Examples to cover the gamut of issues
- Thirty-three page interpretation procedure
- Locus by locus approach
- Pictorial along with calculations and rationale
- Include case specific common scenarios to include victim, suspect and other relevant possibilities
- Highlight and resolve common issues through training



Example 1

All detected alleles are greater than the ST and allelic drop-out was not documented or presumed for this locus. Therefore, all detected alleles are suitable for potential inclusions.

However, to account for the possibility of detected alleles and filtered stutter masking lower level alleles, inconclusive genotypes must still be considered.

Preceding the comparison, a list of inconclusive genotypes can be characterized:

Preceding the comparison, a list of potentially included genotypes can be characterized:

Assume the Victim and Suspect were submitted for comparison. Based on the interpretation of the entire profile, both are included as a donor to this mixture.

Victim: (13, 16) Suspect: (10, 16)


At this locus, the information interpreted as suitable for potential inclusions is consistent with the Victim and the Suspect.

UND = Undetected Allele STT = Stutter

Interpretation of STR DNA Profiles (GlobalFiler®)

 Class:
 Procedure
 Effective Date:
 04/06/2018

 ID:
 Version 5
 Page:
 24 of 33



Example 2

The 12 is less than the ST; therefore, defined as inconclusive.

Based on a resolvability threshold of 4:1, the 8 and the 10 were determined to be unresolved from the 12 allele and therefore, also inconclusive.

$$\frac{446}{127} = 3.26$$

$$\frac{470}{446} = 1.05$$

Therefore, all detected alleles and the resulting genotypes at this locus are inconclusive.


Preceding the comparison, a list of inconclusive genotypes can be characterized:

This locus cannot be used in support of an inclusion or be used to increase the statistical rarity of the profile. The data detected at this locus may only be used to characterize inconclusive genotypes or, if applicable, used for exclusionary determinations.

Interpretation of STR DNA Profiles (GlobalFiler®)

Class:	Procedure	Effective Date.	04/06/2018
ID:	Version 5	Page:	25 of 33





Example 3

The 12, 16, and 19 are less than the ST; therefore, defined as inconclusive.

Based on a resolvability threshold of 4:1, the 14 and the 18 were determined to be unresolved from the 12, 16, and 19 alleles and therefore, also inconclusive.

$$\frac{296}{164} = 1.80$$

$$\frac{385}{296} = 1.30$$

Therefore, all detected alleles and the resulting genotypes at this locus are inconclusive.

Preceding the comparison, a list of inconclusive genotypes can be characterized:

This locus cannot be used in support of an inclusion or be used to increase the statistical rarity of the profile. The data detected at this locus may only be used to characterize inconclusive genotypes or, if applicable, used for exclusionary determinations.


Interpretation of STR DNA Profiles (GlobalFiler*) Procedure Effective Date: 04/06/2018



Version 5

Class:

ID:



Example 5

Page:

27 of 33

All detected alleles are greater than the ST and allelic drop-out was not documented or presumed for this locus. Based on a resolvability threshold of 4:1, all detected alleles are unresolved from one another.

$$\frac{253}{211} = 1.20$$

$$\frac{344}{253} = 1.36$$

$$\frac{428}{344} = 1.24$$

However, based on the number of alleles detected, a minimum of two donors were detected at this locus. Since, the number of donors is less than three, all detected alleles at this locus are suitable for potential inclusions.

Preceding the comparison, a list of inconclusive genotypes can be characterized:

Preceding the comparison, a list of potentially included genotypes can be characterized:

Assume the Victim and Suspect were submitted for comparison. Based on the interpretation of the entire profile, both are included as a donor to this mixture.

Victim: (12, 16) Suspect: (13, 16.3)

At this locus, the information interpreted as suitable for potential inclusions is consistent with the Victim and the Suspect.


Procedures and Training

- Provide a visual
- Use of data and demonstration
- Provide rationale and logic
- Colorful examples to demonstrate gamut of possibilities
- Marry examples to procedures
- Evaluation of procedures though use of a challenge set

Challenge set of 3 complex mixtures

- Three mixtures provided to staff and group of trainees
- Each are 3 person mixtures
- One mixture expected result is use of all loci for statistic, one mixture no loci used for statistic, one mixture with some but not all loci used for statistic
- Permits an evaluation of staff capability
- Evaluates our current protocols
- Provides opportunity for improvement on an individual and group level


Q1 Concatenated Log₁₀ CPI Valuations & Number of Loci Omitted



Mixture 1

- Expected outcome: Mixture profile is too complex for inclusionary purposes
- None of the 40 participants calculated a statistic for this mixture
- A single examiner in training, wrongly interpreted the presence of a major contributor to this mixture
- Provides a teachable moment
- Evaluate training
- Improve understanding
- Produce better, more consistent analysts


Q2 Concatenated Log₁₀ CPI Valuations & Number of Loci Omitted 2



Mixture 2

- Expected outcome: All loci and detected alleles suitable for potential inclusions and statistical estimation purposes (if applicable)
- All 40 participants included and excluded the correct individuals and calculated the stat as expected with an exception of a single participant that omitted one locus
- One examiner, an experienced examiner, chose not to include a locus
- Opportunity to gather information on reasoning and make an appropriate correction
- Notice by the graph there is one "nail" sticking up that needs to be tapped down (or at least evaluated)

Q3 Omitted Loci from CPI - Trainees/Analysts < 6 months



Mixture 3 - Trainees

- Complex mixture containing multiple major contributors suitable for restricted interpretation
- Variation seen among group of trainees, albeit within a generally accepted range
- In the last graph D13 & D2 should have been omitted from being used for potential inclusionary and statistical estimation purposes....D1S should have as well however only 63.6% of the trainees omitted this locus.
- D16 and TH01 were variable (due to potentially low level baseline level activity); either using the locus or omitting it (with appropriate and proper documentation) from the being used for potential inclusionary and statistical estimation purposes was determined to be acceptable
- More documentation needed to describe reason for loci omissions for interpretation and statistic estimation purposes

Take home messages

 Training to include overview, ground truth samples, procedures, good visual case relevant examples, locus by locus demonstration with logic

Use variation of mixtures to evaluate and improve staff

 Continuous improvement means this is a moving target, not a snapshot in time

Be open to change and improvement

Questions?

Thank-you for your kind attention

Ray Wickenheiser

Ray.Wickenheiser@Troopers.NY.gov

(518) 457-1208